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

DISPOSABLE MULTIFUNCTION PLATES ISTRUZIONI FOR USE

PRODUCT DESCRIPTION

A.M.I. ITALIA disposable multifunction pads are made by a couple of pre-gelled self-adhesive electrodes, each equipped with connector for direct connection to the defibrillator/stimulator device.

The rounded shape of the electrodes contributes to limit the risk of current density disuniformity, especially during the defibrillation discharge, so the risk of skin burns. The protection material high insulation together with the safety plug eliminate the unwanted electric micro-shock risks.

A.M.I. ITALIA multifunction pads are available in the following versions: adult, pediatric.

The product is intended to be used by medical staff only in not sterile room.

PACKAGE CONTENTS

A.M.I. ITALIA disposable multifunction plates are packed in a sealed opaque pouch in order to protect gel from light and humidity.

Sales packaging is set by a box containing the pouches and one Instruction for Use sheet.

The product is identified by the selling REF printed on the box.

INDICATIONS

A.M.I. ITALIA disposable multifunction plates are indicated for the following clinical applications:

- external defibrillation
- transthoracic and/or transesophageal synchronized cardioversion
- Temporary transthoracic cardiac stimulation (non-invasive)
- Elettrocardiographic monitoring

A.M.I. ITALIA disposable adult multifunction plates allow the operator to deal efficiently with the treatment of rhythm troubles related with the above mentioned applications, without the risk of any accidental electrocutions related to the normal use of reusable pads.

ATTENTION: it is possible that some of the above indicated clinical applications are not practicable depending on interface system or unit in use.

CONTRAINDICATIONS

- A.M.I. ITALIA disposable adult multifunction plates are not indicated if patient weight is less than 25 kilos or patients are less than 8 years old.
- A.M.I. ITALIA disposable pediatric multifunction plates are not indicated if patient weight is more than 25 kilos.
- Do not apply on irritated or damaged skin.

INSTRUCTIONS FOR USE

For any above described indication, plates must be applied following carefully the instructions of "INSTRUCTIONS FOR APPLICATION" paragraph.

• External Defibrillation, Synchronized Cardioversion: A.M.I. ITALIA disposable multifunction plates are able to transfer to patient the electricity supplied by the defibrillator up to the maximum permitted value 360Joule in the adult version, 130 J in the pediatric version. Indispensable for conversion success, Myocardium critical mass depolarization is possible only if it is crossed by an adequate intensity current: active plate surface is optimized for this aim. For this reason it is opportune to choose the correctly application surface and to apply the plate with the most possible contact surface with the skin. The choice of the energy to be supplied shall be made by the Operator. For pediatric applications the protocol of ERC (European Resuscitation Council) suggests an energy supply of 4J/Kg.

ATTENTION: Every plates couple can support up to 50 defibrillation discharges. Anyway changing them every 24 hours is suggested.

ATTENTION: for transesophageal cardioversion do not exceed 130 Joule for adults and 50J for children.

• Non-invasive transthoracic stimulation: A.M.I. ITALIA disposable adult multifunction plates can be used for non-invasive transthoracic stimulation.

It is appropriate to apply the plates as above described, in order to minimize the stimulation threshold. More over it is necessary to have a good knowledge of the unit components intended for the use and to follow carefully the manufacturer instructions.

ATTENTION: disposable adult multifunction plates must be replaced every 8 hours of stimulation. In case of extended electrostimulation (more than 30 minutes), check if patient skin does not present any irritation sign.

ATTENTION: after an adequate transthoracic stimulation period, it is suggested to keep in consideration to change to an endocavity stimulation.

• ECG monitoring: disposable adult multifunction plates can be used for ECG monitoring

ATTENTION: to continue ECG monitoring during and after the eventual pacing, it is necessary to use a separate ECG electrodes set.

INSTRUCTION FOR APPLICATION

- A.M.I. ITALIA disposable adult multifunction plates can be applied to the patient even in suspicion of a serious arrhythmic pathology.
- Plates positioning points are indicated in figure # 1, 2 and 3 in relation to the intended clinical application. For the best use, please consult the unit manufacturer instruction for use, too.
- Uncover the thorax and prepare the skin. If necessary, remove hairs in excess. Make a light abrasion to the skin to reduce contact impedance.
- Remove eventual residuals using non inflammable detergents. Make sure that the application areas are clean and dry.
- Open the package just immediately before use.
- Gently remove the protective film to expose the adhesive and conductive zones. Do not use the pads if during the detachment of the protective film the product appears damaged (for example the isolating protection of the contact is separated, portions of gel are separated, drops have been formed in the support foam and/or in the pad).

- Apply the plates separately , by pushing the electrode from the top to the bottom with a uniform movement. Plates must be well separated from each other.
- Gently press on the gelled surface and on the adhesive ring to grant a perfect adhesion to the skin and to avoid any air bubble.
- Connect the pads connector to the defibrillator or to the adaptor cable.
- For stimulation on request , connect the electrodes for ECG monitor separately.

SIDE EFFECTS

- Plate adhesive may cause light cutaneous irritation.
- The prolonged transthoracic stimulation or the repeated administration of defibrillation shock may cause more or less noticeable skin reddening according to the supplied energy.
- A lack in adhesion and/or air presence under the electrode may cause burnings.

PRECAUTIONS AND WARNINGS

- Check that the product is compatible to the unit to be used
- Read carefully the instructions for use of the unit to be used
- For pediatrics and for some automatic defibrillators the use of specific energy reduction devices or the use of special shrewdness may be required. A.M.I. ITALIA pediatric electrodes are not equipped with energy reduction systems, therefore please pay particular attention to the energy levels set up on the defibrillator.
- After a prolonged period of transthoracic stimulation, above all with monophasic shaped stimulation waves, the ability of sensing the evoked ECG signal may be reduced. Actually, the consequent increase of the offset continuous voltage of the electrodes can cause the saturation of the signal acquisition and amplification circuits of the connected device. In this case it is necessary to take the evoked signal through a separate set of ECG electrodes.
- Check the expiry date on the plates pack. Do not use after the expiry date.
- Check that packing is complete: do not use in case of damages or visible defects.
- Do not use plates if gel is detached from the support or if it is broken, not homogenous or dry.
- Do not fold the plates and do not put anything on their top.
- Do not use the pads if the connector, or the cable seem to be damaged
- Before discharging energy from defibrillator, be sure that operators are not in contact nor with the patient or any surrounding conductive parts to avoid any accidental electrical shock risk.
- Pay the maximum attention when defibrillators are used in proximity of oxygen source or any other inflammable anaesthetic gas, to avoid any fire or explosion risk.
- Product is not sterile and not sterilizable.
- Product is disposable and to be used for one patient only. Throw it away after use.
- Do not use the plates if they have been removed from the package for more than 24 hours. Plates have to be used within 30 minutes from the protective film removal.

POTENTIAL COMPLICATIONS

- Any complications referred to the use of electrodes aren't expected.

ATTENTION: Defibrillator discharge may cause function irregularity of an implanted pacemaker. After the defibrillator discharge, check if pacemaker is correctly working.

ATTENTION: If the chosen energy level is insufficient the success of the therapy may be compromised. On the contrary, higher levels than necessary could modify the enzymatic order without any evidence of real myocardial damage.

STORAGE

Product must be preserved in its original wrapping in local characterized by environmental conditions, temperature and humidity, specified in the label located on the package. The overlap of weights on the package could damage the product.

PRODUCT LIFE

See the expiry date on sales packaging and on the pouch.

DISPOSAL

Refuses deriving from health structures must be disposed in according to the regulation in force.



WARRANTY AND LIMITATIONS

A.M.I. ITALIA s.r.l. guarantees that the product complies with Directive 2007/47/CE. No responsibility may be ascribed to the producer who shall not be held liable for medical costs or any direct or indirect damage due to lacking function or malfunction of the above product, when used differently from the instruction for use. We recommend to report opportunely any malfunction or defect of the product to A.M.I. ITALIA Quality Assurance Service.

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

Picture 1

- Emergency defibrillation
- Emergency cardioversion of ventricular arrhythmia


Picture 2

- Emergency defibrillation
- Emergency cardioversion of ventricular arrhythmia
- Stimulation

Picture 3

- Cardioversion of atrial arrhythmia



Notes

Electrodes front-front placing (Picture 1) supplies a Lead II trace and this type of placing can be used for cardiac stimulation, defibrillation and monitoring.
The front-rear placing (Pictures 2 and 3) is not recommended for completely automatic and semiautomatic defibrillators (AED)

	European Conformity.		Warning: read the enclosed documentation		Do not reuse
	Production date		Batch number		Use before
	Catalogue Number		Latex free		Keep away from sun light
RH x% y%	Humidity limits		Keep away from humidity		Fragile: handle carefully
	Pcs. per box/pack		Manufactured by		Temperature limits

Note: It is possible that some of these symbols are not included in these Instructions for use.

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